



CM1200B

12-channel Electrocardiograph



Hidden handle




Alphanumeric keyboard;
Back light for button



Memory of 300 ECG records;
USB for data transmission;
PC-ECG management system

Clear professional report



12X1 + Average template + Automatic analysis

✓ 5.7 inch TFT color screen, easy for observation

✓ Internal aluminum alloy rack, shockproof, Anti-interference design



✓ High quality printing system, print is clear and smooth

✓ Internal modular design, easy for maintenance

CM1200B

12-channel Electrocardiograph


COMEN

Share With The World.


Features

- 5.7 inch TFT color screen
- Lead status indication
- Simple & clear keyboard layout, back light support
- Internal memory of 300 ECG records
- High accurate digital filter (EMG filter, AC filter, Drift filter, lowpass filter)
- 122 kinds of diagnosis report
- 120 seconds ECG waveforms review and print
- Support 210mm rolling, Z-folded paper
- USB port for external printer connection & data transmission
- PC-ECG management system
- Working mode: Auto, Manual, Rhythm

Interface



Various printing mode



Printing mode: 12X1, 6X2+1R, 3X4+1R, auto analysis



CM1200B

12-channel Electrocardiograph

Shenzhen Comen Medical Instruments Co., Ltd.

Add: Floor 7, Block 5, 4th Industrial Park of Nanyou, Nanshan District, Shenzhen 518052, China

Tel: +86-755-2640 8879 2641 9446 Fax: +86-755-2643 1232

Website: www.szcomen.com E-mail: info@szcomen.com

CE 0434

