



10 liter
Oxygen Concentrator
introduction





Screen display content

1. Oxygen concentration monitoring function
2. Oxygen flow control
3. Timing function
4. Cumulative time of use



FUSE 5A



power switch


Universal wheel

Light and portable, easy to move. Four universal wheels at the bottom, easy to push.



High efficiency oil - free compressor, strong power, can be 24 hours continuous operation.

Oilless compressor



Molecular sieve

it can effectively control the high level oxygen density

Four layer filtration

- 01 Spongy layer
- 02 Filter box
- 03 Molecular sieve
- 04 Humidification cup





Double filtration makes oxygen fresher.

Details



01

Placement of
oxygen generator
accessories



03

Replace filter
cotton




05



02

Replace filter
cotton



04

Replace filter
cotton



Replace filter
cotton

Technical parameters:

Model number: SD-10W

Max.Flow: 1OL/min

Oxygen density: 93 \pm 3% Running

Noise:s58dB(A)

Dimension: 280*450*555mm

- RATED POWER: 650W
- N.W: 29kgs
- Continuous Working Time: 24 hours

Features:

- Adopts American PSA technology, use Physical method to separate pure oxygen from air;
- Advanced oil-free compressor, save 10% power energy;
- 24 hours continuous working available, 10,000 hours working time warranty;
- 8 years experience and technology make sure the stable and reliable quality;

Useful:

- To strengthen your heart, reducing the risk of heart attacks
- To detoxify your blood, Raise energy level
- To improve strength & build endurance
- To revitalize the skin
- To prevent disease, relieve hangovers
- To help fight off depression

